



PIN	CONNECTION
1	EARTH
2	SOURCE PLATE
3	SEM
4	FILAMENT 1
5	EXTRACTOR
6	SUPPRESSOR
7	RF1
8	REPELLER
9	N/C
10	FILAMENT 2
11	RF2
12	COLLECTOR



PIN	CONNECTION
1	EARTH
2	SOURCE PLATE
3	N/C
4	FILAMENT 1
5	EXTRACTOR
6	N/C
7	RF1
8	REPELLER
9	N/C
10	FILAMENT 2
11	RF2
12	COLLECTOR

ANALYSER VERSION	DIM A	DIM B	DIM C	DIM D
4" SINGLE FILTER	151.70	138.15	156.65	146.70
5" SINGLE FILTER	176.70	163.15	181.65	171.70
DOUBLE FILTER	176.70	163.15	181.65	171.70
TRIPLE FILTER	225.70	212.15	230.65	220.70
TRIPLE FILTER/CHANNELTRON	301.50	287.95	306.45	296.50